

## Non-Esterified Fatty Acids (NEFA) Assay Kit

### Product Information

<b>Cat#</b>	Kit-007
<b>Description</b>	Assay method: ACS-ACOD method
<b>Storage</b>	2–10 °C
<b>Synonyms</b>	NEFA-HA Test
<b>Applications</b>	For the in vitro quantitative determination of free fatty acids in serum. Clinically used as an aid in the diagnosis of endocrine disorders, liver diseases, and diabetes.
<b>Product Overview</b>	<p><b>Kit Components:</b></p> <p>NEFA-HA Test: 2 × 19 mL</p> <p>Free fatty acid calibrator: solution: 4 × 10 mL</p> <p><b>Product Features:</b></p> <p>Ready-to-use format, convenient and efficient</p> <p>Strong method specificity</p> <p>Minimal interference from ascorbic acid and hemolysis; bilirubin may cause slight negative deviation</p>
<b>Notes</b>	Expiration date: 12 months