

## N-Acetyl-beta-D-glucosaminidase (NAG) Assay Kit

## **Product Information**

| Cat#             | Kit-008   |
|------------------|---|
| Description      | Assay method: 4HP-NAG substrate method  |
| Storage          | 2–10 °C   |
| Synonyms         | L-Type NAG  |
| Applications     | For the in vitro quantitative determination of NAG in urine. Clinically used as an aid in the diagnosis of kidney diseases. |
| Product Overview | Kit Components:   |
|                  | L-Type NAG buffer: 2 × 60 mL  |
|                  | L-Type NAG substrate solution: 2 × 20 mL  |
|                  | NAG calibrator: 4 × 3 mL  |
|                  | Product Features:   |
|                  | Most specific indicator for renal tubular function detection  |
|                  | Ready-to-use reagent with stable performance  |
|                  | Morning urine samples not required  |

Fax:1-631-938-8127 45-1 Ramsey Road, Shirley, NY 11967, USA