

N-Acetyl-beta-D-glucosaminidase (NAG) Assay Kit

Product Information

Cat#	Kit-008
Description	Assay method: 4HP-NAG substrate method
Storage	2–10 °C
Synonyms	L-Type NAG
Applications	For the in vitro quantitative determination of NAG in urine. Clinically used as an aid in the diagnosis of kidney diseases.

Product Overview	Kit Components: L-Type NAG buffer: 2 × 60 mL L-Type NAG substrate solution: 2 × 20 mL NAG calibrator: 4 × 3 mL Product Features: Most specific indicator for renal tubular function detection Ready-to-use reagent with stable performance Morning urine samples not required
-------------------------	--