



**Creative Enzymes**

*Diagnostic Enzymes*

## Urokinase

### Product Information

|                         |   |
|-------------------------|---|
| <b>Cat#</b>             | TRA-053   |
| <b>Specification</b>    | 100KIU,500IU,2500KIU  |
| <b>Description</b>      | Urokinase is used in fibrinolytic system research and in vitro thrombolytic reagent development.                    |
| <b>Applications</b>     | Urokinase can be used for the plasminogen activation test and also for the cultivation of three-dimensional fibrin. |
| <b>CAS No.</b>          | 9039-53-6   |
| <b>Synonyms</b>         | Urokinase, from Human Urine   |
| <b>Form</b>             | White or off-white lyophilized powder   |
| <b>Activity</b>         | ≥ 50000 IU/mg   |
| <b>Molecular Weight</b> | 54 kDa/33 kDa   |
| <b>Purity</b>           | ≥ 90%   |
| <b>Endotoxin Level</b>  | < 1.0 EU  |
| <b>Storage</b>          | Store at -20°C for 1 year. Protein Buffer A may also be stored at 4°C.  |

Tel: 1-631-562-8517 1-516-512-3133

Email:[info@creative-enzymes.com](mailto:info@creative-enzymes.com)

Fax:1-631-938-8127

45-1 Ramsey Road, Shirley, NY 11967, USA