



Creative Enzymes

Diagnostic Enzymes

lipoprotein lipase

Product Information

Cat#	EXWM-3460
Similar	LPL
Description	Hydrolyses triacylglycerols in chylomicrons and low-density lipoproteins. Also hydrolyses diacylglycerol.
Form	Liquid or lyophilized powder
CAS No.	9004-02-8
Storage	Store it at +4 oC for short term. For long term storage, store it at -20 oC ~ -80 oC.
Synonyms	clearing factor lipase; diglyceride lipase; diacylglycerol lipase; postheparin esterase; diglyceride lipase; postheparin lipase; diacylglycerol hydrolase; lipemia-clearing factor
Enzyme Commission Number	EC 3.1.1.34
Reaction	triacylglycerol + H ₂ O = diacylglycerol + a carboxylate
Notes	This item requires custom production and lead time is between 5-9 weeks. We can custom produce according to your specifications.

Tel: 1-631-562-8517 1-516-512-3133

Email:info@creative-enzymes.com

Fax:1-631-938-8127

45-1 Ramsey Road, Shirley, NY 11967, USA