

## Taq Antibody (Glycerol-free, 40 U/μL)

## **Product Information**

Cat#	DIA-450
Description	Binds Taq DNA Polymerase reversibly to block polymerase activity at low temperatures preventing non-specific template amplification and primer dimer formation. Denatures during initial DNA denaturation, forming a hot-start Taq enzyme suitable for high-specificity PCR, complex background amplification, low-copy gene detection, high GC content (>60%), and large-scale genomic amplification. Glycerol-free, suitable for lyophilized formulations.
Unit Definition	One unit (U) is defined as the amount of antibody required to inhibit more than 95% of the activity of 1 U Taq DNA Polymerase at 55 °C for 30 minutes.
Storage	–20 °C
Purity	≥ 95%
Contaminants	Free of RNase and endo-/exonuclease contamination
Applications	Compatible with conventional PCR conditions using Taq DNA Polymerase.
Product Overview	Product Features: High affinity for Taq polymerase, effectively preventing nonspecific PCR amplification Rapid hot-start activation: Taq activity fully activated at 95 °C within 30 s
Package	25 μL/125 μL/10 mL

Fax:1-631-938-8127 45-1 Ramsey Road, Shirley, NY 11967, USA