

## AdvSTART Reverse Transcriptase (Glycerol-free)

### Product Information

<b>Cat#</b>	ENZD-MDLyo01
<b>Specification</b>	10000 U; 40000 U
<b>Unit Concentration</b>	200 U/μl
<b>Description</b>	AdvSTART Reverse Transcriptase (Glycerol-free) is a formulation of AccurSTART Reverse Transcriptase developed without glycerol. This modification preserves the exceptional performance and stability characteristics of the original enzyme while enhancing its suitability for lyophilization through the complete removal of glycerol. This product is provided without excipients, allowing for their addition based on specific application requirements.
<b>Notes</b>	Glycerol-free warm-start RT enzyme.
<b>Features</b>	Freeze-drying; Excellent stability
<b>Applications</b>	RT-qPCR is amenable to lyophilization, enabling the development of a simplified one-step assay.
<b>Components</b>	5 × AdvSTART Buffer; AdvSTART Reverse Transcriptase (Glycerol-free) (200 U/μl)
<b>Shipping and Storage Conditions</b>	Store at -30 ~ -15°C, transport at ≤0°C.