



Creative Enzymes

Diagnostic Enzymes

AdvStart II U+ One Step RT-qPCR Probe Kit (for Fast)

Product Information

Cat#	ENZD-RT10
Specification	200 rxns 1000 rxns 5000 rxns
Applications	Diagnostic Pathogen Detection (Respiratory, Intestinal, etc.); Animal Pathogen Detection (African Swine Fever, etc.); Food Testing, etc. (Microbial, etc.)
Description	Suitable for rapid RT-qPCR detection of RNA templates, with fast amplification speed and high efficiency. Contains a dUTP/UDG anti-contamination system, supports multiple amplifications, and ensures detection accuracy.
Features	<ol style="list-style-type: none">1. High sensitivity: Superior to common products on the market, detectable about 1 CT value earlier.2. Rapid amplification: Stable detection within 22 minutes.3. Direct amplification capability: Compatible with a large volume input of nasopharyngeal swab lysis solution 1/2.
Quality Control	High efficiency and sensitivity 22 min fast amplification Compatible with 1/2 volume of lysate
Components	5 × One Step U* Mix One Step U* Enzyme Mix 50 × ROX Reference Dye 1 50 × ROX Reference Dye 2
Shipping and Storage Conditions	Store at -30 ~ -15°C, transport at ≤0°C.

Tel: 1-631-562-8517 1-516-512-3133

Email:info@creative-enzymes.com

Fax:1-631-938-8127

45-1 Ramsey Road, Shirley, NY 11967, USA