



Creative Enzymes

Diagnostic Enzymes

High II One Step RT-PCR Kit

Product Information

Cat#	ENZD-RT21
Specification	50 rxns (50 µl/rxn)
Applications	Virus Detection; Gene Expression Analysis; Endpoint PCR Detection
Description	Endpoint RT-PCR detection for RNA templates, reverse transcription and PCR completed in one tube, improving throughput and reducing contamination risk. Supports up to 10 kb amplification fragments, in 2× Master Mix form, easy to operate.
Features	<ol style="list-style-type: none">1. Reverse transcription and PCR are performed in the same tube.2. Detection sensitivity can reach 1 pg of total RNA.3. Can amplify fragments longer than 10 kb.
Components	RNase-free ddH ₂ O 2 x One Step Mix One Step Enzyme Mix 10 x DNA loading buffer
Shipping and Storage Conditions	Store at -30 ~ -15°C, transport at ≤0°C.

Tel: 1-631-562-8517 1-516-512-3133

Email:info@creative-enzymes.com

Fax:1-631-938-8127

45-1 Ramsey Road, Shirley, NY 11967, USA