



Creative Enzymes

Diagnostic Enzymes

Fluorescent LAMP/RT-LAMP Kit

Product Information

Cat#	ENZD-RT09
Specification	100 rxns 500 rxns
Applications	LAMP/RT-LAMP Isothermal Amplification
Description	Provides 2× fluorescent LAMP/RT-LAMP Master Mix, supports real-time fluorescent detection of DNA or RNA templates. Integrates Bst II Pro enzyme and RTv reverse transcriptase, possessing both high specificity and stability for room temperature preparation.
Quality Control	Excellent LAMP & RT-LAMP amplification performance. High template compatibility. Supports system establishment at room temperature.
Components	2 × Fluorescent LAMP/RT-LAMP Master Mix 50 × LAMP Fluorescent Dye
Shipping and Storage Conditions	Store at -30 ~ -15°C, transport at ≤0°C

Tel: 1-631-562-8517 1-516-512-3133

Email: info@creative-enzymes.com

Fax: 1-631-938-8127

45-1 Ramsey Road, Shirley, NY 11967, USA